

BDSI

www.ukdentalsupplies.com

#whdentalwerk



video.wh.com



Show casing W&H

Special offers
September - December 2024



New
Lisa Mini
Speeds you up!

Restoration & Prosthetics

Contra-angle Synea Power Edition

For the hardest materials



Launch Price
55% Off

The Synea Power Edition is especially designed for heavy-duty applications, such as the removal of decayed materials, cavities, fillings, crowns, and bridges. They are the optimal solution for advanced ceramic material like **zirconia** that is increasingly being used in restorative dentistry.



Synea Fusion WG-900 LT / WG-900 A

- › Contra-angle 1:4
- › Optic / Without light
- › 3x Spray
- › Head Ø 10.7 mm

£ 650

RRP £1,445

£ 585

RRP £1,295



Synea Vision WK-900 LT

- › Contra-angle 1:4
- › Optic
- › 3x Spray
- › Head Ø 10.7 mm

£ 849

RRP £1,699

synea



Contra-angles

High torque and precision

1:5 Speed increasing contra-angles, with FG turbine burs, are recommended for cavity and crown preparation. Due to the maximum speed of 200,000 rpm, aerosol creation is reduced.

1:1 contra-angles are recommended for finishing and polishing.



Synea Fusion WG-99 LT / WG-99 A

- › Contra-angle 1:5
- › Optic / Without light
- › 4x Spray
- › Head Ø 9.5 mm

£ 710

RRP £1,420

£ 635

RRP £1,269



Synea Fusion WG-56 LT / WG-56 A

- › Contra-angle 1:1
- › Optic / Without light
- › Single Spray
- › Head Ø 9.5 mm

£ 562

RRP £1,123

£ 448

RRP £896



Alegria WE-56

- › Contra-angle 1:1
- › Without light
- › Optional External Spray
- › Head Ø 10.0 mm

£ 125

RRP £345



Alegria HE-43

- › Straight Handpiece 1:1
- › Without light
- › Optional External Spray

£ 120

RRP £300



Turbines

High speed, low weight

All W&H turbines have a hygienic head¹ anti-retraction system to prevent cross infection.

¹Exclusive W&H system patented in Europe and the USA

Head Ø

Sprays

No light

Light

Ring LED+

LED+ light



Alegria TE-97/98

- › Without light
- › 3x Spray
- › Head Ø 10.4 / Ø 12.2 mm
- › 16 / 18 W
- › 390,000 / 330,000 rpm
- › Roto Quick® Connection

£ 239

RRP £551



Synea Fusion TG-97 / TG-97 L

- › Without light / LED+
- › 4x Spray
- › MINI Head Ø 10 mm
- › 18 W/390,000 rpm
- › Roto Quick® Connection

TG-97
TG-98 £ 488

RRP £976

TG-97 L
TG-98 L £ 562

RRP £1,123

The TK-94 L microturbine is perfect for hard to access areas and the mouths of children.



Synea Vision TK-94 L

- › LED+
- › 5x Spray
- › MICRO Head Ø 9 mm
- › 14 W/410,000 rpm
- › Roto Quick® Connection

£ 690

RRP £1,458



Air motors & Couplings



Couplings
Available from

£ 101

RRP £203



Air motors
Available from

£ 229

RRP £434

Special Offer Price
When you buy
3x TG 97/98
£ 1,300 Non-optic
£ 1,500 Optic

Our turbines are available to fit other manufacturers couplings.

Prophylaxis & Periodontology

Proxeo Ultra

Scaling can be so easy!

The Proxeo Ultra is available with the unique **Q-Link quick** connection system. In addition, **handpieces are available to fit other manufacturers systems**. The tip changers with integrated torque control prevent overtightening and make changing the tip simple and safe.



GERMAN DESIGN AWARD WINNER 2020



Proxeo Ultra PB-520

- › Compatibility with pacemakers
- › LED ring for ultra-bright illumination
- › Infinitely adjustable power setting
- › Large tank (750 ml)
- › 3 tips (2x universal & 1x Perio)

£ 1,175

RRP £2,205

Compatible with pacemakers

Tested by Graz University of Technology

Full range of tips

The enhanced W&H piezo scaler tip range includes different tips for prophylaxis and periodontology, endodontics, implantology and restoration & prosthetics. The W&H range of tips were developed in close co-operation with end users. The range therefore fulfils the exact requirements of dental practice and implantology.

Tips from

£ 54

RRP £98

Contact W&H for details



W&H Q-Link

W&H / EMS / Mectron

W&H / NSK / Satelec

Proxeo Air Scalers

Regular removal of calculus and plaque leads to a much greater sense of well-being and self-confidence amongst your patients.



Proxeo ZA-55

- › Without power regulation
- › Without light
- › Roto Quick® Connection
- › Tips 1AU, 2AU, 3AU, 1AP

£ 499

RRP £1,249



Proxeo ZE-55 RM/BC

- › Without power regulation
- › Without light
- › Borden/Midwest Connection
- › 2 Tips 1AU

£ 341

RRP £716

Polishing by rotation

Enjoy professional tooth cleaning and give your patients a beaming smile.



Proxeo WP-64 MU

- › Prophy contra-angle 4:1
- › Without light
- › For universal screw-in and snap on (adaptor) cups and brushes

£ 230

RRP £561




Proxeo HP-44 M

- › Handpiece 4:1
- › System Doriot
- › For disposable contra-angles with plastic shank

£ 199

RRP £318



Oral Surgery & Implantology

Surgical Handpieces

High-quality stainless steel with a **scratch-resistant coating** makes the contra-angle handpieces robust and gives them a particularly long service life.



Self-generating LED Surgical Handpieces available.



Contra-angle 20:1 WI-75 E/KM

- › Without light
- › External single spray
- › Not dismantlable

£ 631

RRP £1,262



Contra-angle 20:1 WS-75

- › Without light
- › External single spray
- › Dismantlable

£ 779

RRP £1,559



Contra-angle 20:1 WS-75 L

- › Mini LED+
- › External single spray
- › Dismantlable

£ 852

RRP £1,705



Contra-angle WS-91 / 1:2.7

- › Without light
- › 45° head with 3x spray
- › Dismantlable

£ 734

RRP £1,468



Contra-angle WS-91 L / 1:2.7

- › Mini LED+
- › 45° head with 3x spray
- › Dismantlable

£ 812

RRP £1,625



Straight handpiece S-11 / 1:1

- › Without light
- › External single spray
- › Dismantlable

£ 629

RRP £1,259



Straight handpiece S-11 L / 1:1

- › Mini LED+
- › External single spray
- › Dismantlable

£ 886

RRP £1,772



Handpiece 20:1 Zygoma SZ-75

- › Without light
- › External single spray
- › Dismantlable

£ 1,190

RRP £2,527



Oral Surgery & Implantology

Implantmed Plus SI-1023


Workflow driven solution

Flexibility to meet your needs!



W&H offer a synergistic solution for oral surgery guiding the clinician through the cycle from implementation to evaluation.

You can enjoy the benefits of utilising the modular Implantmed surgical motor and Piezomed piezosurgery unit alongside the Osstell Connect consultation tool and the Osstell Beacon, which ensures predictable assessment of ideal loading.



Implantmed Plus SI-1023 Kit 01 Non-optic/Cable

- Non-optic motor
- Cable Foot Control
- WS-75 20:1 Contra-angle

£ 3,582

RRP £7,163



Implantmed Plus SI-1023 Kit 04 Optic/Wireless

- Optic motor
- Wireless Foot Control
- WS-75 L 20:1 Contra-angle

£ 4,475

RRP £8,950



Wireless foot control for controlling multiple W&H devices

Straight & Contra-angle Handpieces



Osstell Beacon Measure and evaluate the stability of the implant




Implantmed Classic SI-923 Non-optic/Cable

- Non-optic motor
- Cable Foot Control
- WI-75 E/KM 20:1 Contra-angle

£ 3,107

RRP £6,214

Display menu and documentation in your local language



Motors EM-19 EM-19 LC With and without light



Piezomed Classic & Module

"A real practice builder"



Piezomed Classic Optic/Cable

- Piezomed Classic
- Handpiece with LED ring
- Cable Foot Control
- FREE BONE SET Worth £719**

£ 5,633

RRP £9,389



Piezomed Plus Module

- Power 24W. Handpiece with light SA-40 L
- Automatic instrument detection
- Instrument changer and instrument library
- FREE EXTRACTION SET INCLUDED**

£ 3,295

RRP £6,020

Piezomed Instruments

The W&H Piezomed instrument range includes different instruments for bone surgery, periodontology and extraction, for retrograde endodontics and for lateral sinus lift.



Automatic instrument detection



Bone surgery

£ 719

RRP £1,308



Extraction

£ 228

RRP £415



Periodontology

£ 318

RRP £578



Retrograde endodontics

£ 565

RRP £1,027



Lateral sinus lift

£ 516

RRP £938



Preparation of the implant site and crestal sinus lift

£ 1,013

RRP £1,841

Osstell Beacon

Measurement of implant stability



OSSTELL A QUERO COMPANY



Claim **FREE SmartPegs***



Osstell Beacon

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and avoid rejection in high-risk patients.

£ 1,989

RRP £5,123

Osstell SmartPegs

The SmartPeg measures the resonance frequency when it is attached to an implant.

£ 83

RRP £113



Osstell Connect FREE App for Osstell Beacon

Osstell Connect is designed to document quantified data which ensures traceability, whilst complementing data in your Patient Management System.

*Warranty only for Implantmed Plus SI-1023 kits.

Sterilization, Hygiene & Maintenance

PremiumCare



Gives you the flexibility to select a plan that best suits the individual needs of your practice.

The standard PremiumCare package offers a flexible plan allowing you to add optional extras as you please or the **PremiumCare Plus or Flex**** plan, which include fully comprehensive service packages and the required pressure vessel inspection and validation testing at no extra cost.

PremiumCare

Maintenance plan for all W&H sterilizer models, with a maximum limit of 4,000 cycles over the 3 year term*.

PremiumCare Plus & Flex**

A fully comprehensive maintenance plan for all W&H sterilizer models, including PVI certification, unlimited cycles and validation in line with manufacturer's or current HTM 01-05 guidance (model dependent).

Our Technical Service Team are here to support you nationwide!



Optional
Not included
Included

Services	Premium Care	Premium Care Plus & Flex**
All routine service visits for three years	✓	✓
All service parts and scheduled sundries	✓	✓
All breakdown calls	✓	✓
All breakdown parts including labour	✓	✓
All engineer travel in the UK Excluding Highlands, Islands and ferries etc,	✓	✓
Free on-site operator training	✓	✓
Free technical software upgrades	✓	✓
Unlimited Cycles - for the duration of the plan	⊘	✓
PVI Certification - in line with Pressure Systems Safety Regulations 2000:128	⊘	✓
Periodic Test Logbook	⊘	✓
Validation in line with manufacturer's or current HTM 01-05 guidance (model dependent)	✗	✓
Installation and Commissioning in line with current guidelines	✗	✓
Helix Test Kit***	✗	✓
Extra Discount	✗	✓



Internal cleaning & lubrication

Increase your handpieces lifetime and reduce the costs of repairs!



quick connect
by W&H

Assistina Twin

- Two process chambers allowing loading of a 2nd handpiece while the 1st is processed
- Superior Maintenance & Lubrication
- Intuitive and easy to use
- A flexible range of adaptors available

£ 1,400
RRP £3,108



quick connect
by W&H

Assistina One

- Internal cleaning & lubrication of one handpiece in only 18 seconds, with Assistina One the aerosols stay in the HEPA filter thereby releasing clean air.

£ 1,020
RRP £2,105



Assistina One includes Care set and Quick RM and ISO adaptors



Service Oil F1
Ask your dealer about W&H Handpiece Oil's



3x EliReel Included

Seal2 Starter Pack
Universal Pouch Sealing
12mm hermetic seal
3 x EliReel Rolls FREE

£ 656
RRP £1,208

W&H EliReel sealing pouches

From / per roll
£ 13
RRP £21



Osmo
Supplies multiple sterilizers simultaneously with top quality demineralized water.
Automatic water filling
Most cost-effective solution
Increase sterilizer service life
Up to 6 devices simultaneously

£ 899
RRP £1,498



Seal Check
Pack of 100 sealing check

£ 29
RRP £60



Lisa Safe Label Printer & QR Barcode Reader
Quick and easy tool for linking the sterilization cycle and dental products used for treatment

Barcode reader
£ 142
RRP £227

Lisa Safe printer
£ 499
RRP £777

*An additional cycle tariff will be levied if the 4,000 cycle limit is exceeded within the 3 year period ** PremiumCare Flex Plans are only available for certain models *** Only applies to B Type Sterilizers

Sterilization, Hygiene & Maintenance

Cleaning & Disinfection

The new Teon washer disinfecter is easy to use and offers an exceptional standard of cleaning, disinfection and reprocessing, maximising efficiency and increasing the effectiveness of the decontamination protocol.



BePro Cleaner TE
Tri-enzymatic non-foaming detergent for automatic cleaning in washer disinfecter.

£ 44
RRP £54



BePro Rinse
For rapid, stain-free drying during automated reprocessing of dental instruments.

£ 14
RRP £17

Teon thermal washer disinfecter
 › Validated cleaning of all instruments
 › Traceability and connectivity through Ethernet module
 › 1 Year Warranty

Teon with Standard Starter Set
£ 6,800
RRP £9,071

Teon with Handpiece Injector Starter Set
£ 8,100
RRP £10,257

BePro WDD Starter Pack
 › BePro WDD Cleaner TE – 5000 ml
 › BePro WDD Rinse – 1000 ml
 › Canister opener

£ 69

134°C
SSS

Sterilization at its Best!

With a wide range of sterilizers, W&H is able to offer the perfect option to meet individual needs.



Lyla S type Sterilizer
New concept sterilizer
 › Mobile App available via upgrade for remote data storage
 › Future-proof upgradability to a type B sterilizer using the new Activation Code concept

Unique upgrade flexibility with **Activation Codes**
 › B Cycle
 › Remote data storage
 › Traceability
 › All in One

£ 3,318
17 l. RRP £4,608

£ 3,696
22 l. RRP £4,847

All W&H Sterilizers offer:

- › Intuitive user interface and touchscreen display for a smoother user experience
- › Cycle traceability available thanks to Wi-fi connection and built in data logger
- › Ergonomic design, with modular feet for flexible positioning
- › Easy to clean surfaces and accessible water › Automatic Water Fill function
- › Pre-programmable vacuum and helix (for B Type Sterilizers) test function to run automatically out of surgery hours



Lara B type Sterilizer
 › Mobile App for Real Time Remote Monitoring
 › Eco Dry technology automatically adapts the drying time to the mass of the load
 › Unique optional upgradability to increase: Performance and Traceability or to add a Fast Cycle

Unique upgrade flexibility with **Activation Codes**
 › Performance**
 › Traceability
 › Fast Cycle
 › All in One
 ** Only for Lara 17 and 22 litres

£ 4,559
17 l. RRP £6,642

£ 5,121
22 l. RRP £7,461

Lara XL
£ 5,648
28 l. RRP £8,230



Lisa B type sterilizer
 › Full extended traceability thanks to EliTrace
 › Mobile App for Real Time Remote Monitoring
 › Eco Dry+ technology automatically adapts the drying time to the mass of the load
 › Elisense user interface generates status, temperature and other cycle information
 › Cycle times from 13 minutes for specific loads
 › State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use

£ 5,467
17 l. RRP £8,205

£ 6,178
22 l. RRP £9,286

Lisa MINI
£ 2,870
3 l. RRP £4,100

£ 3,200
5 l. RRP £4,600

NEW Lisa MINI B type sterilizer
 › Cycle times from 10 minutes
 › Integrated vacuum pump
 › Small and easily stackable: W 44cm x H 20cm x D 67cm
 › Built in water tank and automatic water filling - no external storage
 › For up to 22 instruments, offers complete traceability

*Based on average 2kg load for immediate use

**with fast cycle upgrade

BePro

Professional Cleaning and Disinfection for All Surfaces



Soaking bath

Soaking bath for the immersion of dental instruments in the decontamination phase.

£ 37
RRP £44

Disinfectant R 1l.

Ready-to-use solution for the disinfection and cleaning of rotary dental instruments.

£ 13
RRP £15

Disinfectant I 2l.

Powerful disinfectant concentrate for medical and dental instruments.

£ 38
RRP £46

Disinfectant Spray 1l.

Surface disinfectant for the rapid disinfection and cleaning of non-invasive medical devices.

£ 13
RRP £15



Disinfectant Wipes

120 wipes for the disinfection of non-invasive medical devices.

£ 9 £ 7
Tube RRP £10 Refill RRP £8

Meet the Sales Team



Jon Bryant

National Sales and Marketing Manager

07917 886326

jon.bryant@wh.com



Richard Mansbridge

Area: London, South East, Anglia, East Midlands and Humberside

07775 792628

richard.mansbridge@wh.com



Caroline Caine

Area: West Midlands, Manchester, Merseyside South West England, Isle of Man, Channel Islands, Wales, Northern Ireland and Eire

07786 515331

caroline.caine@wh.com



David Barrass

Area: Scotland, Humberside, North of England and North West

07795 315089

david.barrass@wh.com

W&H (UK) Ltd

6 Stroud Wood Business Centre,
Park Street, St Albans, AL2 2NJ

t +44 (0) 1727 874990

f +44 (0) 1727 874628

office.uk@wh.com

www.wh.com

To place an order
or for more information
contact BDSI

01226380400

Offers valid from 1st September to 31st December 2024. Information is correct at time of publication and relates to W&H products and services available in the UK and Ireland. Prices are quoted in £ sterling and exclude VAT.

Savings are based on RRP and rounded to the nearest £.

Free goods are not transferrable and no cash equivalent is available.

Terms and Conditions apply. E. & O. E.