



Technical specifications

Product data

Camera Specifications

Sensor	CMOS
Resolution	1920 x 1080px
Distortion	-8%
Frame per second (FPS)	25 - 30
Aspect ratio	16/9
Lighting	8 white LEDs
Focus adjustment	Fixed focus
Image capture	Capture button or footswitch (optional)
Field Of View Diagonal (FOVD)	70°
Cable length	3m
USB output	2
Consumption	2.5W
Handpiece dimensions	L: 200,9; W: 26.9; H: 23.3mm
Handpiece weight	35g
International Protection Marking (IP)	Dust & waterproofness protection IP53

Configurations

Computer requirements

Specs

CPU	Intel® Core i5 4 core 2GHz Gen 8 or greater and Apple M1 or greater
RAM	8GB
USB version	2.0 or greater
USB chipset	Intel or NEC®/RENESAS®
Operating system	from Windows 10 PRO 64 bits or macOS® Ventura

